

# BLOPAD™

## Caring for Wounds

**BIOPAD is an equine collagen primary wound dressing used to treat minor bleeding, open ulcers, and skin lesions.**

- **100% pure native equine type I collagen**

**BIOPAD has No fillers**, just pure collagen to stimulate proliferation of fibroblasts and for formation of granulation tissue.

- **Highest collagen content on the market**

With 5 times the standard amount of collagen, BIOPAD does not breakdown as quickly as the competitors' thin dressing, sending more collagen to the wound site.

- **Easy to use**

BIOPAD can be used in a hyperbaric chamber, with NPWT and can be used to fill large voids.

- **Cuts costs when used with SilverStream®**

When supplemented with SilverStream, BIOPAD eliminates the need to use a higher priced silver/collagen dressing.

- **Soft gel dressing**

BIOPAD transforms into a soft gel that allows contact within the entire wound bed.



# BIOPAD™

Caring for Wounds

**BIOPAD**  
is also  
**Reimbursable**

## Instructions for Use:

- Before BioPad application, properly prepare wound bed per standard protocols removing necrotic tissue where present.
- Apply BioPad directly to the wound using sufficient amount to completely cover the wound bed.
- For dry wounds, hydrate the BioPad dressing with sterile saline or SilverStream® to initiate the gelling process.
- Secure the BioPad with a non-adherent secondary dressing of choice to maintain moisture.
- The gel that is formed by interaction with the wound exudate is biodegradable and does not need to be removed at each dressing change, unless it is necessary to cleanse the wound for some other reason (e.g. presence of infection).

**For more Information  
contact us at:**

**1-800-726-2308**

## Presentation and size:

BioPad is packaged in a transparent waterproof blister pack and is sterilized by gamma irradiation

2" x 2" - 3 per box

**Order code:** 132622 **HCPCS code:** A6021



ANGELINI

### CONTACT US

**Angelini Pharma Inc.**

8322 Helgerman Court - Gaithersburg - MD 20877 -USA

1 (800) 726-2308

[www.Angelini-us.com](http://www.Angelini-us.com)

20-011-01 Rev. 01 2/14